



The Hudson
range

Operation/Care & Warranty Manual



PRIDE[®]
MOBILITY PRODUCTS LTD.
LIVE YOUR BEST[®]

Fine Quality British Built Reclining Chairs and Complementary Lounge Sofas



Thank you for your purchase

We are very proud of our brand. Our Chairs and Sofas represent the very best of British Built Motion Furniture.

The Hudson range was established in 2016 with a commitment to produce the finest quality Chairs and Sofas with the very best mechanisms and framework. Each piece of furniture is thoughtfully designed, traditionally handmade and then subjected to intensive quality inspection by skilled craftsmen here at our factory in Nottinghamshire.

The Hudson range of furniture reflects our passion for quality. We only source the finest materials to produce long lasting sumptuous and comfortable seating and our frames are made from selected sustainable hardwoods and are built to last.

We believe that our Chairs and Sofas are not merely functional objects, they are part of our wellbeing and what make our home life comfortable and warming; it gives us great pleasure helping you sit comfortably and to be able to get in and out of practical uplifting furniture.

This manual is to help you look after your furniture. If used in accordance with these guidelines, your furniture could last a lifetime.



Operation/ Care & Warranty Manual



Assembly

All Recliners / lift recliners come in two pieces for easy installation purposes.

Recliner/Lift recliners & reclining sofas


Simply slide the brackets on the backs of the chair over the male brackets situated at the back of the seating area, a gentle click noise means the chair back has locked into place (Fig 1).

If a swivel locking bracket is used, simply follow the same instructions but instead of an automatic click lock a manual swivel lock is needed; simply push the swivel arms into place on the male brackets to lock (Fig 2).

To take apart, simply reverse the operation by releasing the spring clips (click lock) and pulling up the back from under the wings one lock at a time-or pull out the swivel arms (swivel lock) and pull the back up from under the wings.



• (Fig 1) - Click Lock



• (Fig 2) - Swivel Lock

Luxury Lift & Recline Chairs - Single Motor Operated



- Returns chair from recline and lifts the chair
- Reclines chair
- Standby button - Turns on the handset

Safety feature

Please note: our handsets can be disconnected by unclipping the linking ports inside the pocket of your chair.

Luxury Lift & Recline Chairs - Dual Motor Operated




Standard Dual Motor Handset

- Returns foot plate, lifts the chair
- Reclines foot plate
- Returns backrest
- Reclines backrest
- Standby button turns on the handset

Standard controller Dual Motor - VivaLift! Mechanism

Control Operation



Standard controller with powered headrest and lumbar support chair - VivaLift! Mechanism

These lift chairs feature a 10" powered headrest, powered lumbar support, dual motor design and USB charging functionality.



Cleaning & Appearance

To keep your furniture looking as new as possible please follow these important instructions:

- Avoid long periods of sunlight exposure
- Lightly vacuum regularly
- Upholstery must be cleaned using products advised by professionals. Do not use potentially harmful liquids or chemicals

Appearance of your furniture

This chair or sofa is built with the highest quality fillings and foams and is designed for longer periods of sitting compared to standard domestic non motion furniture. Therefore it is inevitable that your chair/sofa will start to mould to your shape or sitting position. Please note that seating quilts on reclining furniture will require brushing from time to time using your hands and 'dressing' after long periods of sitting or reclining. Gently brush out any creases or friction lines, and regularly plump up the back cushion after long periods of sitting to avoid similar distortion and to keep the appearance as new as possible.

Battery Back Up System

All reclining furniture (electric) is supplied with a battery back up transformer. Two 9 volt batteries are fitted in the transformer which will provide reassurance that in the unlikely event of a power failure, the chair will return to a sitting position if stuck in an elevated or reclined position (batteries will need replacing once emergency system has been used).

WARRANTY

Your warranty period specified by your retailer, covers the following;

- Electrical and mechanical parts, (mechanisms, motors, handsets, transformers, plugs, etc)
- Wooden framework
- Foam build - including seat foams
- Seating springs

All the above are covered providing the product remains with the original purchaser and the product has been used for domestic purposes. It is important that this product has been used in accordance with these operating and care guidelines and has not been misused or damaged.

If your recliner fails to operate then please follow these instructions:

- Check all leads are connected from the chair to the transformer and from transformer to mains electric plug (light should be illuminated on the transformer)
- Check that no cables have been trapped or damaged by the chair's operation,

Once all the above have been checked, and if the chair is still not working, please call your retailer.

Please note: Fabric or leather is not covered by the full warranty period, this comes with a standard 12 month guarantee.

Do's & Don'ts

Things that can and will affect warranty claims upon inspection.

Don'ts

- ✗ Sit on the arms of any motion/static furniture
- ✗ Recline your furniture into any object situated close by
- ✗ Attempt to modify any of the furniture's specification
- ✗ Attempt to repair or replace any electrical malfunction
- ✗ Store any solid objects under bed mechanisms

Do's

- ✓ Keep your motion furniture away from any object that may obstruct its movement
- ✓ Report an electrical fault or any other service issue with your retailer as soon as possible
- ✓ Make sure children are not left alone with any motion furniture product
- ✓ Wipe away any spillages on upholstery immediately with a damp cloth
- ✓ Treat your furniture with general care - it could last a lifetime!

Warranty Claims

Any warranty claims need to be given to your retailer at the earliest possible moment, warranty claims will be dealt with upon inspection at our factory and need to have an incident/service report with the product along with the original purchase order, invoice date and also the SP serial number (situated on the base mechanism at the rear of the chair or bed).



The Hudson
range

Your Notes

Your Notes



The Hudson range



32 Wedgwood Road, Bicester, Oxon, OX26 4UL
www.pride-mobility.co.uk Email: sales@pride-mobility.co.uk
Tel: +44 (0) 1869 324600 Fax: +44 (0) 1869 323070



@pridemobilityuk

Ver: 1.3
March 2019

PRIDE
MOBILITY PRODUCTS LTD.
LIVE YOUR BEST*

Fine Quality British Built Reclining Chairs and Complementary Lounge Sofas