



LION

OUTDOOR ROLLATOR

LION ROLLATOR









FRAME COLOUR



Anthracite



The Mobilex Lion is designed for use outdoors where the ground surface is varied. It is perfect on off road trails, over cobblestones and gravel pathways where its large soft wheels can ensure a comfortable ride. Although robust, it is still relatively lightweight, and offers a high degree of stability and ease of movement.

The push handles are height-adjustable with numerical markings and a simple lock. It features a textile seat, padded back strap, cane holder and a unique locking system that stops it reopening once folded.

The large 9" smooth-running rear wheels are equipped with brakes for optimum speed and movement control, and are also lockable for parking. The mesh bag has a load capacity of 5kgs and can be removed and used as a convenient shopping bag.

Specifications	ML53	ML62
Overall Width	68cm	68cm
Overall Length	76cm	76cm
Product Weight	9.4kgs	9.7kgs
Width Folded	26cm	26cm
Seat Height	53cm	62cm
Seat Width	56cm	56cm
Height of Push Handles	77cm – 90cm	90cm – 112cm
Max. User Weight	200kgs	200kgs
Wheels Front Rear	285mm x 44mm 235mm x 44mm	285mm x 44mm 235mm x 44mm
Max. Bag Load	5kgs	5kgs

OPTIONAL ACCESSORIES



WATERPROOF BAG

PADDED SEAT





ONE HANDED BRAKE WHEELS

Sold and maintained by: