

clinell®

Alcohol Wipes

Surface disinfection
of non-invasive
medical devices.



Rapid, high quality alcohol wipes

Ultra-quick drying, pharmaceutical grade 70% isopropyl alcohol wipes.

Disinfection wipes, which leave no smear, for surface disinfection of non-invasive medical devices.

Proven efficacy



High-quality durable wipes containing pharmaceutical grade 70% isopropyl alcohol. Our wipes are proven to have bactericidal and virucidal action.

Unique packaging



Wipes are dispensed one at a time, so the top of the next wipe is not dry, which commonly occurs with tub dispensers. This eliminates the infection control risk of dry contaminated wipes.

Less plastic. Less waste



Our flow wrap packaging is environmentally friendly as there is much less packaging to dispose of compared to standard alcohol tubs.

Superior quality



Clinell Alcohol Wipes contain 25gsm non-woven spun bond, which is wetter & thicker than competing alcohol wipes allowing for larger area disinfection with fewer wipes, saving time and money.

For the surface disinfection of non-invasive medical devices:

- Ultrasound probes
- Pulse oximeter probes
- Doppler probes
- Stethoscopes
- ECG machines
- Blood pressure machines
- Blood pressure cuffs
- Radiology equipment
- Incubators
- Monitors

Antimicrobial efficacy



Proven to kill at least 99.999% of microorganisms.

Including MRSA, *E. coli*, *Pseudomonas*, *Enterococcus*, *Candida* and *Aspergillus* according to EN1276, EN13727, EN14561, EN13624 & EN14562.

	Microorganism example	EN test
Bacteria	<i>Acinetobacter baumannii</i>	EN13727
	<i>Enterococcus hirae</i>	EN13727
		EN 16615
	<i>Escherichia coli</i> K12	EN13727
	Methicillin resistant <i>Staphylococcus aureus</i>	EN13727
	<i>Pseudomonas aeruginosa</i>	EN13727
		EN16615
	<i>Staphylococcus typhimurium</i>	EN13727
Fungi	<i>Staphylococcus aureus</i>	EN13727
		EN16615
Mycobacteria	<i>Aspergillus brasiliensis</i>	EN13624
	<i>Mycobacterium avium</i> (2 year aged)	EN14348
	<i>Mycobacterium terrae</i> (2 year aged)	EN14348
Viruses	Norovirus	EN14476
	Murine coronavirus	EN14476
	Vaccinia virus	EN14476
Yeast		EN13624
	<i>Candida albicans</i>	EN16615

Ordering information

Product	Unit of issue	Order code	NHSSC code (UK only)
Alcohol Wipes Large 150	Pack of 150	CAW150L	VJT839

For more information, speak to your GAMA Healthcare Area Manager or visit www.clinell.com



CE 0050 Class IIa Medical Device
Conforms to the Medical Device Directive, 93/42/EEC.

Use disinfectants safely. Always read the label and product information before use. Always follow the medical device manufacturer's cleaning and disinfection guidelines prior to use with Clinell Alcohol Wipes.

GAMA Healthcare Ltd.,
The Maylands Building, Maylands Avenue,
Hemel Hempstead, HP2 7TG, UK.

T: +44 (0)20 7993 0030
E: info@gamahealthcare.com
www.gamahealthcare.com

JBN220359