

Don't leave home without it

**FreeRider**<sup>®</sup>  
working together

**FR168-3**  
**KNIGHTSBRIDGE 3**



**NHC**  
Mobility

**300 lbs**

**4 mph**



# FR168-3

## KNIGHTSBRIDGE 3



### SPECIFICATIONS:

Overall length	115 cm/45 inches
Overall width	61.5 cm/24.2 inches
Ground clearance	9.5 cm/4 inches
Maximum load	136 kgs/300 lbs
Range up to	32 km/20 ml
Total weight (excl batteries)	49.53 kgs/109 lbs
Heaviest Part	23.5 kgs/40lbs
Wall to wall turning radius	96 cm/38 inches
Maximum speed	6.4 km/4 mph
Battery	36 Amp

#### Electronics

State of art, fully proportional speed control for safe operations.

#### Drive Systems

Rear Wheel direct drive via transaxle

#### Motor

24 volt DC permanent magnet, totally enclosed for outdoor use.

#### Brakes

Automatic dynamic regenerating system with magnetic solenoid parking brake.

#### Free wheel facility

#### Automatic battery charging systems

### FEATURES:

- Distinctive styling
- Easy to use controls
- Comfortable seat rotates, slides, height adjustable
- Lift-up arms for easy access
- Sturdy Shopping basket
- Infinitely adjustable steering column
- Maintenance-free sealed batteries
- Front and rear lights
- Styled 10 inches pneumatic wheels front
- Styled 11 inches pneumatic wheels rear
- Solid tyres (option)



**CONTACT YOUR DEALER TODAY**

Our policy is one of continual improvement.

We reserve the right to change specification without prior notice.

All figures stated are approximate and do not constitute accurate specifications.

The colours shown are for general guidance and are as close as can be achieved by modern printing process.