

Royale 3 & 4

ROYALE



The Drive Royale 3 and 4 are two of the finest mobility scooters available today, both models have an impressive performance; powered by 70Ah batteries they have an excellent maximum range of 32 miles. These stylish scooters also feature a modern illuminated LCD display, providing Speedo, Battery Gauge, Clock, Odometer, Temperature, Reverse Light and Intelligent Diagnostics.

- › Top speed of 8mph*
 - › Maximum range of 32 miles on full battery charge*
 - › Stylish modern design
 - › Comfortable rotating seat with sliding and reclining facility
 - › Height adjustable headrest to suit user requirements
 - › Armrests are angle and width adjustable and flip up for easy side transfers
 - › Automatic speed reduction when steering for added safety
 - › All round suspension ensures a smooth ride and optimal handling
 - › Tiller adjustment by lever to suit user requirements and assist easy transportation
 - › Digital dash display with user friendly controls
 - › High level indicators and lighting to improve outdoor visibility
- › at dusk or when dark
 - › Twin rear view mirrors for extra visibility
 - › Large 40cm (16") alloy wheels with black pneumatic tyres for a smoother ride
 - › T-Bar style handlebars
 - › Anti roll back system safety device prevents scooter from rolling backwards on hills
 - › Freewheel facility allows movement of scooter without switching the motor on
 - › Rear anti tip wheels as standard
 - › Lap belt included for added safety
 - › Optional rear pannier box
 - › Available in Red, White and Black
 - › Recommended MK batteries: 2 x 12V 79Ah
 - › Weight capacity 185kg (29st)
- * Subject to conditions of use

Colour Options

- Red
- White
- Black

“The result is an exceptional driving experience”

CODE	DESCRIPTION	LENGTH	HEIGHT	WIDTH	MAXIMUM GRADIENT	MAXIMUM WEIGHT CAPACITY
HW001	Royale 3 Wheel	158cm (62")	135cm (53")	73cm (29")	10°	185kg (29st)
HW002	Royale 4 Wheel	158cm (62")	135cm (53")	73cm (29")	10°	185kg (29st)

