

Airium

Motion
healthcare

HEAVIEST PART JUST 13KG!



Lithium-Ion
Battery Technology



Lightweight Aluminium
Frame



25 Mile Range



25 Stone Weight Limit



AVAILABLE IN:



Plum



Teal

OTHER FEATURES



SOFT GRIP DELTA
HANDLEBAR



COMFY ROTATING
SEAT



LARGE WHEELS WITH
SUPER SOFT TYRES



REAR SUSPENSION

LENGTH	112cm / 44"
WIDTH	62cm / 24.4"
HEIGHT	93.5cm / 36.8"
FOLDED LENGTH	N/A
FOLDED WIDTH	N/A
FOLDED HEIGHT	N/A
SEAT WIDTH	46cm / 18"
FRONT WHEEL	Solid 22.5cm / 9"
REAR WHEEL	Solid 22.5cm / 9"
TOTAL WEIGHT	104.3lbs / 47.31kg
HEAVIEST PART	28.6lbs / 13kg
BATTERY WEIGHT	15.4lbs / 7kg
MAX USER WEIGHT	25stone / 160kg
RANGE	25miles / 40km
SPEED	4mph / 6.4kph
BATTERY	25.2v 31.5ah Lithium
MOTOR	24v 250w
GROUND CLEARANCE	5.2cm / 2"
MAX GRADIENT	8°
TURNING RADIUS	145cm / 57.1"
BOX LENGTH	130cm / 51"
BOX WIDTH	55cm / 21"
BOX HEIGHT	53cm / 20"
PACKAGED WEIGHT	108.3lbs / 49.12kg



VISIT www.motionhealthcare.co.uk

CALL 0844 257 8150