



Oak



Walnut

The **Alerta Gas Lift Overbed Table** features an easy lift handle to adjust the table height with minimal effort. It is finished in oak or walnut and comes with brakes on the two rear castors. The table has a laminate top and comes with a wheelchair accessible base design.

TECHNICAL SPECIFICATION

- > Wheelchair access base
- > Braking castors
- > Max weight capacity: 15kg
- > Height adjustment range: 80-105cm
- > Table top dimensions: 80 x 40cm

PRODUCT CODES

ALT-GL-O - Alerta Gas Lift Overbed Table, Oak

ALT-GL-W - Alerta Gas Lift Overbed Table, Walnut

KEY FEATURES AND BENEFITS



Rear castors



Easy lift height adjustment



Won't topple over even with 15kg on top!



Oak without castors



Walnut with castors

The **Alerta Overbed Table** provides easy access to everyday conveniences for users seated in a bed, chair or wheelchair. It is finished in oak or walnut and comes with or without castors. The table has a laminate top with two adjusting points which allows both height and tilt adjustment and comes with raised edges to prevent items slipping off.

TECHNICAL SPECIFICATION

- > Max weight capacity: 10kg
- > Height adjustment range: 64-97cm
- > Table top dimensions: 61 x 40cm
- > Max Tilt angle: 34°

PRODUCT CODES

ALT-751O - Alerta Oak Overbed Table, No Castors

ALT-751W - Alerta Walnut Overbed Table, No Castors

ALT-751OC - Alerta Oak Overbed Table, With Castors

ALT-751WC - Alerta Walnut Overbed Table, With Castors

KEY FEATURES AND BENEFITS



Castor option



Height & tilt adjustment



Won't topple over even with 10kg on top!