

Titan LTE Powerchair

Compact,
manoeuvrable
and easily
transportable



Red & blue
interchangeable
colour panels



The Titan LTE is a compact, lightweight and highly transportable chair which is simple and easy to use due to its user friendly programmable controller. This powerchair splits into 3 sections for easy handling and storage when not in use.

www.drivedevilbiss.co.uk

drive
 **DeVilbiss**
HEALTHCARE

Combining performance and portability with a powerchair that is nimble, responsive and comfortable. It is simple to operate and disassemble for transportation.

Interchangeable Colour Panels



Red & Blue Included
Code: TITAN201

Features

- ▶ Maximum range of 8 miles on a full battery charge*
- ▶ New & easy to use Dynamics LiNX controller system
- ▶ 2 x 250W Motors
- ▶ Designed for indoor and outdoor use
- ▶ Compact & manoeuvrable for use in tight spaces & busy environments
- ▶ Splits easily into 3 pieces for easy storage
- ▶ Convenient carry handle on the rear base unit
- ▶ Detachable & width adjustable armrests to suit comfort of user
- ▶ Seatbelt included for additional safety and security
- ▶ Height adjustable seat allows chair to be easily adjusted to suit comfort preferences
- ▶ Interchangeable colour panels
- ▶ Swivel seat facilitates easy transfers on and off the chair
- ▶ Freewheel Mode to allow use of the powerchair without switching the motor on
- ▶ Padded upholstery for additional comfort and support
- ▶ 20cm (8") front castors and 23cm (9") rear wheels with solid puncture proof tyres provide durability and low maintenance
- ▶ Height adjustable & flip-up footplate to suit users leg length and assist with easy transfers
- ▶ Rear anti-tip wheels as standard
- ▶ Comes with batteries and off board charger

*Subject to conditions of use

Specifications

Width:	60cm (24")
Length:	88cm (34.5")
Height:	87cm (27.5")
Max Speed:	4mph*
Max Climbing Angle:	6°
Ground Clearance:	5cm (2")
Front Wheels:	8"
Rear Wheels:	9"
Batteries:	12V 20 Ah x 2
Total Weight:	54.4kg (44lb)
Weight Capacity:	115kg (18st)



New intelligent Dynamics LiNX controller



Splits into 3 sections for easy handling and storage



Easy to carry detachable components

Leaflet Code: LL414

If you would like to view this product please contact your local stockist.

We have a continuous product improvement policy and consequently reserve the right to amend design and specification without prior notice. All sizes and weights stated are nominal. All errors and omissions excepted. Due to the limitations of photography and the printing process colours shown may not be 100% accurate.