

# HAMPTON ACTIVE MATTRESS



**Product Code:**

HAMHAMP001

**Description:**

- The Hampton is Harvest Healthcare's mid-range Active Replacement Mattress
- 20 cells made of Nylon with a TPU lining to reduce noise when cells inflate and deflate
- Pillow section for patient stability
- All cells are interchangeable for ease of repair and cleaning
- 4-way stretch welded PU cover which is vapour permeable and water resistant and has a Crib 5 fire rating
- Compatible with all of Harvest Healthcare's pumps allowing for a range of maximum patient weights

**Mattress Dimensions:**

2000 x 900 x 200mm  
2200 x 900 x 200mm (Also available on request- Comes with Pump 2)

**Weight:**

8kg

**Max. Patient Weight:**

Harvest 1Pump: 25 Stone (158kg)  
Harvest 2 Pump: 29 Stone (184kg)  
Harvest 3 Pump: 40 Stone (254kg)

**Pressure Relief Risk:**

Very High

**Pump Technical Specification:**

Please see separate pump specification sheet

**Cleaning Guidance:**

Cover should be wiped down with hot water and mild detergent. If heavily soiled, the cover may be machine washed up to 65°C and tumble dried. For disinfection use sodium hypochlorite solution max 1000ppm. Use no bleach, biological or phenolic-based cleaners. Dry completely before storing.

**Warranty**

In terms of Harvest Healthcare Guarantee and Warranties the Harvest Hampton mattress is guaranteed for three years from date of purchase.

**Exclusions:**

1. Removable covers are guaranteed for one year
2. Damage to fabric covers through incorrect use and penetration by sharp instruments is excluded.

**Supplier Contact Details:**

t: +44(0)1709 377172      e: [sales@harvesthealthcare.co.uk](mailto:sales@harvesthealthcare.co.uk)      w: [www.harvesthealthcare.co.uk](http://www.harvesthealthcare.co.uk)  
Sheaf house, Bradmarsh Way, Bradmarsh Business Park, Rotherham UK, S60

**Important Safety Notice**

Harvest Healthcare Ltd supplies equipment and appliances to cover a range of medical conditions. Specific choice of equipment is the sole responsibility of the clinical practitioner or carer overseeing individual patient care. Harvest Healthcare Ltd supplies a range of products such as foam mattresses, active air mattresses, fall protection mats and bed rail protection bumpers. These products are usually covered in vapour permeable PU fabric that could cause risk of suffocation in the event of prolonged facial contact. Mattress risk categories are to be used as a guide only. Risk assessments specific to user must be carried out prior to the use of all products. Products and data continually undergo development and may therefore deviate from the details given.

**Version 3: March 2020**