

Technical Specifications

The Scooterpac Canopy



Frame Material		Powder Coated Aluminium	
Canopy Material	Waterproof Nylon	Cover Material	Waterproof Nylon
Window Material		Perspex	
Maximum Wind Speed		40mph	

	Standard	XL
--	----------	----

Open Height ¹	1200mm	1200mm
Open Width	600mm	600mm
Open Depth	1150mm	1280mm
Folded Size	H830 x W630 x D180mm	H960 x W630 x D180mm
Weight	12kg	13kg

¹Varies depending on scooter manufacturer and model.
The Scooterpac Canopy requires a standard mobility scooter accessory bar to be fitted. If your scooter does not feature the accessory bar an adaptor is available. To check if your mobility scooter is compatible please see www.scooterpac.com/check.



Absolute freedom

Stay warm and dry all year round.
For complete freedom, contact your
Scooterpac dealer to have yours
installed today.

www.scooterpac.com
info@scooterpac.com

+44 (0)330 124 1733


scooterpac®
Lifestyle mobility products

Designed in
Great Britain

