

# RIVIERA

Bathlift

## Riviera Bathlift

*"The Riviera has seen many developments over the years and continues to be popular"*



CE

For many people, at some point in their life, bathing becomes difficult. The Riviera bathlift provides a fresh approach to this problem and enables a person to enjoy bathing again with confidence.

**drive**<sup>TM</sup>  
[www.drivemedical.co.uk](http://www.drivemedical.co.uk)

The Riviera lifts, lowers and reclines easily and quietly and is designed to fit in most baths. It is easy to install, simple to operate and offers a person excellent comfort, safety, reliability and convenience.



Hygienic sealed surface covers

## A brilliant, reliable & safe solution to all your bathing needs...

- > **Hygienic** - The backrest design with smooth, sealed surfaces is more resistant to dirt, making a bath lift easier and faster to clean.
- > **Easier** - The storage battery compartment is located in the centre directly behind the backrest, thus making it a lot easier to insert or remove the battery.
- > **Moving and handling** - The bathlift is easily dismantled into two pieces, both of which include an integral handle for moving and handling.
- > **Practical** - The unlocking mechanism of the drive unit features a red unlock button - one simple push with your thumb is sufficient to fold and dismantle the bath lift.
- > **Safer** - Due to an EMERGENCY-OFF button, which is integrated in the hand control, a more powerful storage battery offering even longer operating times, and a unique safety shutdown function.
- > **Space-saving** - Flat backrest design for greater leg room.
- > **Attractive** - The slim line design makes the bath lift even more attractive and adapts harmoniously to its surroundings.
- > **Reliable** - Lifetime guarantee on frame and two years on battery and handset..

## RIVIERA BATHLIFT

CODE: 430100252 (White)  
430100312 (Blue)

Total Width: 70cm (27½") (seat including side flaps)	Total Length: 86cm (34") (at lowered backrest)	Seat Depth: Approx. 45cm (18")	Seat Width: 37cm (14½")
--	--	--------------------------------------	----------------------------


Seat Height: 8cm - 42cm (3" - 16½")	Total Length of Seat: 58cm (23")	Height of Backrest: 56cm (22")	Width of Backrest: 32cm (12½")
---	-------------------------------------	-----------------------------------	-----------------------------------

Backrest Reclines: 12° - 40°	Width of Baseplate: 31cm (12")	Length of Baseplate: 54.5cm (21½")	Duration of Lifting/Lowering: 25 seconds
---------------------------------	-----------------------------------	---------------------------------------	---

Max Capacity Weight: 140kg (22st)	Drive Unit Weight: 7.5kg (17lb)	Frame Weight: 8kg (18lb)	Total Weight with Storage Battery: 16.5kg (37lb)
--------------------------------------	------------------------------------	-----------------------------	---

Leaflet Code: LL121

If you would like to view this product please contact your local stockist, details below.



Easily dismantled



Easy access battery



Easy to use handset