

UltraFine Straight Leg

The Oxford® UltraFine Straight Leg sling is made from a cushioned and breathable 3D woven 'Spacer' material with high air permeability and pressure redistribution properties. Integrated, low-friction rip-stop nylon to the leg sections ensures easier application.

- Available with or without head support
- Detachable leg and shoulder straps
- Colour coded binding tab on top section to denote size
- Coloured attachment loops promote correct/safe positioning
- Modesty loop helps protect patient dignity during transfer
- Can be customised to fulfil patient specific needs


















Chair to Chair



Bed to Chair/Chair to Bed



Chair to Commode/Toilet



Transfer to Bath



Floor to Bed/Chair (Head Support recommended)



Rehab/Assisted Standing

Size	Part No. (Without HS)	Part No. (With HS)	SWL
XS	SL1561	SL1821	227 kg
S	SL1562	SL1822	227 kg
M	SL1563	SL1823	227 kg
L	SL1564	SL1824	227 kg
XL	SL1565	SL1825	284 kg







