

**DASH**<sup>TM</sup>

# Dash Super HD Wheelchair

The heavy duty chair suitable for up to 46 stone user weight.

UP TO  
**46st**  
USER WEIGHT



**FOLDABLE**

Foldable wheelchair frame to enable easy storage, transport or maintenance.



**WIDE SEAT**

Extra wide seat for increased space and comfort. Available in 26", 28" and 30".



**ADJUSTABLE FOOTRESTS**

The footrests can be height adjusted to suit the user. They can also be removed if required.

Dimensions	26"	28"	30"
Total width	880mm	930mm	980mm
Total length	1250mm	1250mm	1250mm
Total height	980mm	980mm	980mm
Seat width	660mm	710mm	760mm
Seat depth	500mm	500mm	500mm
Backrest height	470mm	470mm	470mm
Folded width	360mm	360mm	360mm
Front wheels	8"	8"	8"
Rear wheels	24"	24"	24"

Other technical data	
Frame	Steel
Front wheel	PU
Rear wheel	PU with aluminium rim
Handgrip	PVC
Seat and Backrest	Nylon
Armpad	PU
Footplate	Steel
Brake	Steel



### Weight info

Chair lifting weight 45kg

Max user weight 295kg

## Key features

- ✔ Up to 46 stone (650lb) user weight
- ✔ Foldable frame
- ✔ Height adjust footrests
- ✔ Steel parking brake
- ✔ Puncture proof tyres

### For enquiries:

T: 0845 146 0600

E: sales@rhealthcare.co.uk

F: 01246 453881

[www.rhealthcare.co.uk](http://www.rhealthcare.co.uk)



RHealthcare

It's ability that matters